Fișa suspiciunii de plagiat / Sheet of plagiarism's suspicion	Indexat la: 0106/07
--	--------------------------------

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
Cage System Surgical Technique

10 Steps

- 1 Preoperative Planning and X- Ray Analysis
- 2 Patient Preparation
- 3 Discectomy and Decompression
- 4 Distraction
- 5 Selection
- 6 Solis Bonegraft Harvesting Set Option
- 7 Bonegraft Preparation
- 8 Overdistraction
- 9 Positioning
- 10 Axial Compression

Cage System Surgical Technique

1 Preoperative Planning and X-Ray Analysis


Conventional x-ray templates are available with the Solis Cervical Cage System to assist in the selection of implant sizes.

Overlay the x-ray template on 15% magnification front and lateral x-ray films to determine the appropriate implant size (x-ray yield in average a 15% magnification).

1 - Properly sized implants should fit within the confines of the anterior cortex, posterior cortex and uncovertebral joints on the front x-ray.

2 - Overlay the templates on the lateral x-ray. Measure the height of the adjacent non-degenerative discs until a cervical cage template that matches is identified.

The identified size gives an estimation of the Solis cage that should be used at the degenerated level.

2 Patient Preparation

Place the patient in a supine position, stabilized in a head holder. The patient's head higher than his feet to reduce the venous pressure. A slightly rotated and extended head position to the other side gives a better surgical approach. Too much rotation though gives too much tension of the sternocleidomastoideus muscle.

Prepare and drape the patient in the usual manner for anterior cervical decompression and fusion surgery. The iliac crest on the same side may be prepared.

Use intraoperative radiographic imaging to identify the affected discs and make an incision that optimizes the exposure of the procedure.

3 Discectomy and decompression


The discectomy is complete when all disc material has been removed.

Care shall be taken to remove all the bone that causes the compression of the nervous structures.

Cage System Surgical Technique


4 Distraction

a With the Pin Driver first insert one of the distraction pins in the cervical vertebrae above or beneath the disc to be operated.

b

Place the Pin Guide on the pin and insert the second distraction pin through the sleeve. The Pin Guide allows the placing of the two pins in a parallel position.


Cage System Surgical Technique

4 Distraction


C Place both arms of the Distractor on the pins.

NB: If required the Solis cage system provides two Distractor Blades.

d

Perform a complete anterior cervical decompression by turning forward the wing nut of the Distractor.

Cage System Surgical Technique

5 Selection

a Connect the Trial Cage to the Implant Holder by turning the handle.

In order to prevent the cage mispositioning during the insertion,


- the Implant Holder features an etching which indicates the cranial position,
- and the Trials as well as the final implant have a lateral hole which corresponds to the ball on the tip of the Implant Holder.

NB: The size of the Trials corresponds to the Solis cage without the spikes.

b

The Trial Cage is then placed in the intra-discal space. Try if sufficient compression arises on the inserted trial. If not try the next larger size trial cage. Consult your preimplantation measurements and compare these results with your operative findings.


Cage System Surgical Technique

6 Solis Bonegraft Harvesting Set Option

a To take bonegraft in a minimal invasive way a Bonegraft Harvesting Set is provided and composed by:


b

Cage System Surgical Technique


7 Bonegraft preparation

- a** After connecting the Solis Cage to the Implant Holder as already done with the trials, use the special Bonegraft Compactor.

b

By just turning the detachable base, the Bonegraft Compactor allows the placing of the two depths of the Solis cage: 12 and 14 mm.

Cage System Surgical Technique


8 Overdistraction

NB : In order to accomodate the spikes of the implant, the disc space should be slightly overdistacted.

9 Positioning


Insert the cage into the disc space in the midline of the spine and release the Implant Holder.

Take a lateral x-ray to ensure the Solis cage is in the right position in the front and lateral view.


10 Axial Compression

Push on the button and turn the wing nut backward to release the Distractor and perform axial compression.

Remove the Solis cage distraction device by just taking out the Distractor and then the pins with the Pin Driver.

Cage System Surgical Technique

1 Implants

6741204	SOLIS Cervical Cage 12/4
6741205	SOLIS Cervical Cage 12/5
6741206	SOLIS Cervical Cage 12/6
6741207	SOLIS Cervical Cage 12/7
6741208	SOLIS Cervical Cage 12/8
6741209	SOLIS Cervical Cage 12/9

6741404	SOLIS Cervical Cage 14/4
6741405	SOLIS Cervical Cage 14/5
6741406	SOLIS Cervical Cage 14/6
6741407	SOLIS Cervical Cage 14/7
6741408	SOLIS Cervical Cage 14/8
6741409	SOLIS Cervical Cage 14/9

2 Trials

8741204	Trial Cage - 12/4
8741205	Trial Cage - 12/5
8741206	Trial Cage - 12/6
8741207	Trial Cage - 12/7
8741208	Trial Cage - 12/8
8741209	Trial Cage - 12/9

8741404	Trial Cage - 14/4
8741405	Trial Cage - 14/5
8741406	Trial Cage - 14/6
8741407	Trial Cage - 14/7
8741408	Trial Cage - 14/8
8741409	Trial Cage - 14/9

3 Instruments

874011	Distractor
874002	Distraction Pin 10mm
874003	Bonegraft Compactor
874004	Implant Holder
874005	Pin Driver
874006	Bonegraft Harvesting Kit
874007	Pins guide
874008	Distractor blade
874009	Distraction pin 15mm

4 Containers

874000	Sterilization Container
874000A	Sterilization Micro-Container

5 Templates

874010	Radiographic Template
--------	-----------------------

stryker
SPINE

EU Operations
Z. I. Maricot - 93 810 Cestas - FRANCE
Phone : +33 (0)5 57 97 06 30
Fax : +33 (0)5 57 97 06 31

US Operations
59 Route 17 South Allendale,
New Jersey 07401 - USA
Phone : +1 201 825-4900
Fax : +1 201 790 8108

CE 0459